

Bioengineering 278: Magnetic Resonance Imaging Laboratory

Winter 2012

Lab 6

1. **Generate an interferogram.** An interferogram (in the context of MRI) is an image in which two or more coherence pathways are present in the image and are allowed to interfere with one another. In this exercise, you will generate an interferogram between a spin echo and an FID (gradient echo) generated by the refocusing pulse. You will visualize the interference by applying a shim offset in the X direction. The shim offset should perturb the phase of the gradient echo but not the spin echo. The shim offset will cause the spin echo and gradient echoes to interfere with one another in a spatially dependent manner. Place the pink phantom in the birdcage coil. Use the pulse sequence spep, which is a spiral pulse sequence that is homegrown. Prescribe a single slice axial spin echo sequence which will use default flip angles (90-180). Set CVs: dda=4, reps=1, rf1type=1, pw_rf1=3200, pw_rf2=3200, cyc_rf1=2, cyc_rf2=2, slwid180=1, nl=16. Scan, and verify that you get an image. This is the reference image.
 - a. Calculate what combination of flip angles (rf1 and rf2) will generate a spin echo and a gradient echo from rf2 of equal amplitude (neglecting T1 decay). The solution may not be unique. (3 points)
 - b. Adjust ia_rf1 and ia_rf2 to achieve these flip angles, and scan. The default crusher gradients around the refocusing pulse should destroy the FID, so you should only see the spin echo. Calculate the expected signal intensity relative to the reference image, scan, and compare your results with the predicted values. (3 points)
 - c. Eliminate the effect of the crusher gradients by setting CV:zcrush=-4. Verify using the oscilloscope that the FID from rf2 should not be crushed. Change the value of the X shim by 40 units, and rescan. From the spacing of the stripes, calculate the gain of the manual X shim adjustment (in G/cm per unit shim offset). (3 points)
 - d. Calculate the expected profile across the stripes and compare with the acquired profile. (6 points)
2. **Resonance offset sensitivity in FIESTA (balanced SSFP, FISP).** Place a uniform phantom in the scanner. Scan the phantom using a single axial slice. Using the FIESTA sequence, determine the gain of the manual X shim adjustment (in G/cm per unit shim offset), and compare with the same measurement from part 1. (5 points)