Bioengineering 208 Magnetic Resonance Imaging


Winter 2007 Lecture 1

E. Wong, BE208, UCSD Winter 2007

Topics

- Review of NMR basics
- Hardware Overview
- Quadrature Detection

Boltzmann Distribution



$$\frac{\text{Number of Spins Down}}{\text{Number of Spins Up}} = e^{-\Delta E/kT}$$

Ratio = 0.999990 at 1.5T!!!

Corresponds to an excess of about 10 up spins per million

E. Wong, BE208, UCSD Winter 2007

Slide Credit: T.T. Liu

Magnetization Vector


Vector sum of the magnetic moments over a volume.

For a sample at equilibrium in a magnetic field, the transverse components of the moments cancel out, so that there is only a longitudinal component.

Equation of motion is the same form as for individual moments.

$$\mathbf{M} = \frac{1}{V} \sum_{\substack{\text{protons}\\ \text{in } V}} \mu_i$$

$$\frac{d\mathbf{M}}{dt} = \gamma \mathbf{M} \times \mathbf{B}$$




http://www.easymeasure.co.uk/principlesmri.aspx

E. Wong, BE208, UCSD Winter 2007

Slide Credit: T.T. Liu

RF Excitation



At equilibrium, net magnetization is parallel to the main magnetic field. How do we tip the magnetization away from equilibrium?

B₁ radiofrequency field tuned to Larmor frequency and applied in transverse (*xy*) plane induces nutation (at Larmor frequency) of magnetization vector as it tips away from the *z*-axis.

- lab frame of reference

http://www.eecs.umich.edu/%7Ednoll/BME516/

E. Wong, BE208, UCSD Winter 2007

Slide Credit: T.T. Liu

Bloch Equation


$$\frac{d\mathbf{M}}{dt} = \mathbf{M} \times \gamma \mathbf{B} - \frac{M_x \mathbf{i} + M_y \mathbf{j}}{T_2} - \frac{(M_z - M_0) \mathbf{k}}{T_1}$$
Precession
$$\omega_0 = \gamma \mathbf{B}$$
Transverse
Relaxation
Transverse
Relaxation
Relaxation

 \mathbf{i} , \mathbf{j} , \mathbf{k} are unit vectors in the x,y,z directions.

Slide Credit: T.T. Liu


Relaxation: Z-component

$$M_z(t) = M_0 + (M_z(0) - M_0)e^{-t/T_1}$$



Relaxation: Transverse Component

$$\begin{split} M &= M_x + jM_y \\ dM/dt &= d/dt \big(M_x + iM_y \big) \\ &= -j \big(\omega_0 + 1/T_2 \big) M \\ M(t) &= M(0) e^{-j\omega_0 t} e^{-t/T_2} \quad \omega_0 = \gamma \mathbf{B} \end{split}$$




E. Wong, BE208, UCSD Winter 2007


Slide Credit: T.T. Liu

Relaxation: Summary

- 1) Longitudinal component recovers exponentially.
- 2) Transverse component precesses and decays exponentially.



Gradient Echo




E. Wong, BE208, UCSD Winter 2007

Slide Credit: T.T. Liu

Small Tip Angle Example

$$\begin{split} B_{1}(t) &= B_{1}rect\left(\frac{t-\tau/2}{\tau}\right) \\ M_{r}(\tau,z) &= jM_{0}\exp(-j\omega(z)\tau)\int_{0}^{\tau}\exp(j\omega(z)s)\omega_{1}rect\left(\frac{s-\tau/2}{\tau}\right)ds \\ &= jM_{0}\exp(-j\omega(z)\tau/2)F_{1D}\left(\omega_{1}rect\left(\frac{t}{\tau}\right)\right)\bigg|_{f=-(\gamma/2\pi)G_{z}z} \\ &= jM_{0}\exp(-j\omega(z)\tau/2)\omega_{1}\tau\sin c\left(f\tau\right) \\ &= jM_{0}\exp(-j\omega(z)\tau/2)\omega_{1}\tau\sin c\left(\frac{\gamma G_{z}\tau}{2\pi}z\right) \end{split}$$



E. Wong, BE208, UCSD Winter 2007

Slide Credit: T.T. Liu

Hardware Overview

Three fields:

- •Main Field (B₀)
 - •Polarize Spins
- •Gradient Fields (G_[XYZ])
 - •Map space into frequency
 - $\overline{\partial [XYZ]}$
- •RF Fields (B₁)
 - •Change the latitude or zenith angle of ('tips') spins

E. Wong, BE208, UCSD Winter 2007


Main Field (B₀)

How do we decide on B₀?

$$\Delta E = \gamma h B_0$$

 $M_0 \propto \Delta E$

... Bigger is better!











3T Human @UCSD.

7T Rodent @UCSD

7T Human @U.Minn. 9.4T Human @UIC


Main Field (B₀)

Energy in a Magnetic Field:

$$E = \frac{1}{2\mu_0} \int B^2 dV$$

For B=3T over $1m^3$:

$$E = \frac{1}{2(1.25 \times 10^{-6})} 9 = 3.6 MJ$$
 = dropping a 1000Kg car from 360m high



During a quench, R goes from 0 to $\sim 100\Omega$, $I\sim100A$, so $P=I^2R\sim1MW$

A quench can boil off 3.6MJ/2.5KJ/l=1400l of Helium in 3.6MJ/1MW ~3.6s!!!

E. Wong, BE208, UCSD Winter 2007

Main Field (B₀)

Wavelength (λ) of RF:

In Vacuum:


5m @ 60MHz (1.5T)


1m @ 300MHz (7T)

In Brain:

12cm @ 300MHz (7T)¹

When λ is not large compared to object, standing waves form. This is referred to as Dielectric Resonance. RF inhomogeneity during receive is fixed by scaling, but RF transmit inhomogeneity is much more difficult to address.






¹Vaughan et al, MRM 46 p24 2001

Gradient Fields

$$G_X = \frac{\partial B_Z}{\partial X}$$
 $G_Y = \frac{\partial B_Z}{\partial Y}$ $G_Z = \frac{\partial B_Z}{\partial Z}$



- •Shortest soft tissue $T_2^* \sim 1$ ms
- •For 0.2mm resolution in 1ms:

$$G = \frac{K_{\text{max}}}{\gamma T} = \frac{(0.5/0.2mm)}{(4257Hz/G)(1ms)} \approx 5G/cm$$

• To fill 1m³ with 5G/cm gradients in 0.2ms requires:

$$P = \frac{E}{T} = \frac{1/2\mu_0 \int B^2 dV}{T} \approx \frac{1/2\mu_0 (B_{RMS} (5G/cm))^2 (1m^3)}{0.2ms} \approx 500KW \longrightarrow \frac{\text{About 3 simultaneous}}{\text{Rolling Stones concerts}}$$

- Modern gradient systems are also up against dB/dt limits for peripheral nerve stimulation (~50T/s)
- For diffusion or ultrashort T₂* imaging, more G would help a lot

E. Wong, BE208, UCSD Winter 2007


RF Fields

How big do RF fields need to be?

- •Shortest soft tissue $T_2^* \sim 1 ms$
- •To flip spins by 90° (0.25 rotations) in 0.2ms:


$$B_1 = \frac{0.25}{\gamma T} = \frac{0.25}{(4257Hz/G)(0.2ms)} \approx 0.24G$$

• RF power absorption by the body is a complex function of frequency, conductivity, and geometry, but at 0.24G, approximately 200W/Kg are deposited in human tissue at 3T. Thus, for a 100Kg person, the RF system must supply 20KW of deposited power, or about 40KW of total power, assuming 50% losses to the coil, cabling, reflections, and radiation.



Quadrature Reception

Original quadrature detection: separate coils and A/D for I and Q




E. Wong, BE208, UCSD Winter 2007

http://www.easymeasure.co.uk/principlesmri.aspx

Quadrature Reception

Because ω_0 >>bandwidth, 2 coils are not needed for phase detection. However, second coil does increase SNR




E. Wong, BE208, UCSD Winter 2007

http://www.easymeasure.co.uk/principlesmri.aspx

Quadrature Reception

Digital quadrature requires 2x faster sampling, but eliminates I/Q imbalance, and what happens to a DC offset in the A/D?



Quadrature Reception

Summary:

- 1. Quadrature RF coil is NOT needed to detect MR signal phase.
- 2. Quadrature RF coil improves SNR by $\sqrt{2}$
- 3. Digital quadrature detection:
 - 1. Eliminates I/Q imbalance
 - 2. Moves DC offset in ADC to edge of image
 - 3. Requires 2x higher sampling rate than separate I/Q